

Oil Red O (for frozen sections), Kit

Reference: AP0368



1 of 2

INTENDED USE AND PRESENTATION:

For *in vitro* diagnostic use.

AP0368. 100 Tests.

APPLICATIONS:

The Oil Red O (for frozen sections) kit was designed to demonstration of lipids on 5 µm frozen sections.

Fixation: it is recommended the use of Saline formalin or Baker's fixative in order to maintain the phospholipids.

The interpretation of the stain results is the full responsibility of the user. Any experimental result must be confirmed by a medically established diagnostic product or procedure.

PRODUCT COMPOSITION:

The Oil Red O (for frozen sections) kit is composed by:

Reagent	Name	Volume	Storage
A	Absolute propylene glycol	30 mL	15-25 °C
B	Oil red O solution	30 mL	15-25 °C
C	Propylene glycol 85% solution	30 mL	15-25 °C
D	Mayer's hematoxylin solution	30 mL	15-25 °C

METHODS AND PROCEDURE:

Principles of the method: Oil Red O solution stains selectively cellular lipids while hematoxylin solution stains nuclei.

Specimen: Frozen sections.

Procedure time: 25 minutes.

Staining procedure:

- 1) Rinse the slides containing the frozen sections in distilled water.
- 2) Place in **Reagent A** for 2 minutes.
- 3) Stain in **Reagent B** for 16 minutes.
- 4) Differentiated in **Reagent C** for 1 minute.
- 5) Rinse in two changes of distilled water.
- 6) Stain in **Reagent D** for 15 to 60 seconds .
- 7) Rinse throughly in several changes of distilled water.
- 8) Mount with mounting medio.

See our web site at www.gennova-europe.com for detailed protocols ancillary reagents and support products.

EXPECTED RESULTS:

During the reaction of the different solutions of the kit with the tissue sections, resulting in the following stains:

Lipids	Red
Nuclei	Blue

REQUIRED MATERIALS BUT NOT SUPPLIED:

All reagents, materials, and laboratory equipment for this procedure are not provided with this kit. This includes adhesive slides and cover slips, positive and negative control tissues, xylene or adequate substitute, ethanol, distilled H₂O, pipettes, Coplin jars, glass jars, moist chamber, histological baths, mounting materials, and microscope.

Buffered solutions for and other auxiliary reagents are available from Gennova Scientific.

STORAGE AND STABILITY:

Store at 15-25 °C until the expiration date printed on product label. Do not use after the expiration date. Keep the containers tightly closed. After the first opening, the product is usable until the expiry date stated on label of intact product. Waste from solutions must be disposed under the procedure of hazardous substance. If the product is stored under different conditions from those stipulated in these technical indications, the new conditions must be verified by the user.

Gennova Scientific guarantees that the product will maintain all of the described characteristics from the production date until the expiration date, as long as the product is stored and used as recommended. No other guarantees are provided. Under no circumstances is Gennova Scientific obliged to cover damages caused by use of this reagent.

TROUBLESHOOTING:

If unusual staining is observed or any other deviations from the expected results, please read these instructions carefully, and if this



Catalog number



Batch code



In Vitro diagnostic medical device



Temperature limitation



Expiration date



Test number



Manufacturer



See instruction for use



Gennova Scientific, S.L.
C/ Johann Gutenberg, 4F. Pol. Ind.
El Cádizamo I • 41300 San José
de La Rinconada • Sevilla, SPAIN
Teléfono: +34 954 150767
Fax: +34 955 266494

info@gennovalab.com
www.gennova-europe.com

Oil Red O (for frozen sections), Kit

Reference: AP0368



2 of 2

does not solve the problem, please contact Genova Scientific's technical support department or your local distributor.

PRECAUTIONS:

Use only by qualified personnel. Read with attention the information written on the label (dangerous symbols, risks and safety phrases). Consult always the safety data sheet (MSDS) where the information about the risks of the preparation, precautionary measures during use, first aid and disposal are available. Use proper protective equipment in order to avoid contact with reagents and samples in the eyes, skin, and mucosal tissues. MSDS is available upon request.

PERFORMANCE CHARACTERISTICS:

Genova Scientific has performed studies to evaluate the functioning of this product, concluding that the product is both specific and sensitive for recommended use.

BIBLIOGRAPHY:

Pearse AGE. Histochemistry Theoretical and Applied, 2 nd ed., London.

Bancroft JD, Gamble M: Theory and Practice of Histological Techniques, Churchill Livingstone Elsevier, sixth ed., Philadelphia

F01IT04_V1R0212_AP0368_English



Catalog number



Batch code



In Vitro diagnostic medical device



Temperature limitation



Expiration date



Test number



Manufacturer



See instruction for use



Genova Scientific, S.L.
C/ Johann Gutenberg, 4F. Pol. Ind.
El Caamo I • 41300 San Jose
de La Rinconada • Sevilla, SPAIN
Telefono: +34 954 150767
Fax: +34 955 266494

info@genovalab.com
www.genova-europe.com