

Weigert-Van Gieson quick method **(for elastic fibres and connective), Kit**

Reference: AP0384



1 of 2

INTENDED USE AND PRESENTATION:

For *in vitro* diagnostic use.

AP0384. 100 Tests.

APPLICATIONS:

The Weigert-Van Gieson quick method (for elastic fibres and connective) kit was designed to show elastic fibers, connective tissue and nuclei at the same time in tissue sections. Recommended for vascular pathology.

The interpretation of the stain results is the full responsibility of the user. Any experimental result must be confirmed by a medically established diagnostic product or procedure.

PRODUCT COMPOSITION:

The Weigert-Van Gieson quick method (for elastic fibres and connective) kit is composed by:

Reagent	Name	Volume	Storage
A	Lithium carmine solution	30 mL	15-25 °C
B	Weigert's resorcin fuchsin solution	30 mL	15-25 °C
C	Van Gieson picrofuchsin solution	30 mL	15-25 °C

METHODS AND PROCEDURE:

Principles of the method: The method is based on the affinity towards elastic fibers displayed by the dye resulting from a reaction between resorcin and basic fuchsin in the presence of ferric chloride. Since this is not absolutely specific method, other structures such as collagen and basal membranes might be stained. Therefore it is essential to differentiate carefully in order to obtain a selective marked staining of elastic fibers.

Specimen: Paraffin-embedded tissue samples should be used.

Procedure time: 60 minutes.

Staining procedure:

- 1) Deparaffinize and hydrate for paraffin section.
- 2) Put the slide in distilled water.
- 3) Put on the section 2-5 drops of **Reagent A** for 10 minutes.
- 4) Wash quickly in running tap water.
- 5) Put on the section 2-5 drops of **Reagent B** for 40 minutes.
- 6) Differentiate in alcohol 95°.
- 7) Put on the section 2-5 drops of **Reagent C** for 5 minutes.
- 8) Rinse in distilled water.
- 9) Dehydrate in alcohols of increasing clear in xylene and mount preparation.

See our web site at www.gennova-europe.com for detailed protocols ancillary reagents and support products.

EXPECTED RESULTS:

During the reaction of the different solutions of the kit with the tissue sections, resulting in the following stains:

Elastic fibers	Black
Background	Grey uncolored
Cytoplasm, erythrocytes, muscle fibers	Yellow
Collagen fibers	Red purple
Nuclei	Red

REQUIRED MATERIALS BUT NOT SUPPLIED:

All reagents, materials, and laboratory equipment for this procedure are not provided with this kit. This includes adhesive slides and cover slips, positive and negative control tissues, xylene or adequate substitute, ethanol, distilled H₂O, pipettes, Coplin jars, glass jars, moist chamber, histological baths, mounting materials, and microscope.

Buffered solutions for and other auxiliary reagents are available from Genova Scientific.

STORAGE AND STABILITY:

Store at 15-25 °C until the expiration date printed on product label. Do not use after the expiration date. Keep the containers tightly closed. After the first opening, the product is usable until the expiry date stated on label of intact product. Waste from solutions must be disposed under the procedure of hazardous substance. If the product is stored under different conditions from those stipulated in these technical indications, the new conditions must be verified by the user.



Catalog number



Batch code



In Vitro diagnostic medical device



Temperature limitation



Expiration date



Test number



Manufacturer



See instruction for use



Genova Scientific, S.L.
C/ Johann Gutenberg, 4F. Pol. Ind.
El Caamo | • 41300 San Jose
de La Rinconada • Sevilla, SPAIN
Telefono: +34 954 150767
Fax: +34 955 266494

info@gennovalab.com
www.gennova-europe.com

**Weigert-Van Gieson quick method
(for elastic fibres and connective), Kit**

Reference: AP0384



2 of 2

Genova Scientific guarantees that the product will maintain all of the described characteristics from the production date until the expiration date, as long as the product is stored and used as recommended. No other guarantees are provided. Under no circumstances is Genova Scientific obliged to cover damages caused by use of this reagent.

TROUBLESHOOTING:

If unusual staining is observed or any other deviations from the expected results, please read these instructions carefully, and if this does not solve the problem, please contact Genova Scientific's technical support department or your local distributor.

PRECAUTIONS:

Use only by qualified personnel. Read with attention the information written on the label (dangerous symbols, risks and safety phrases). Consult always the safety data sheet (MSDS) where the information about the risks of the preparation, precautionary measures during use, first aid and disposal are available. Use proper protective equipment in order to avoid contact with reagents and samples in the eyes, skin, and mucosal tissues. MSDS is available upon request.

PERFORMANCE CHARACTERISTICS:

Genova Scientific has performed studies to evaluate the functioning of this product, concluding that the product is both specific and sensitive for recommended use.

BIBLIOGRAPHY:

Humberstone, G.C.W. e Humberstone, F. D.: An elastic tissue stain. J. Med. Lab. Technol., 26, 99, 1969
Moore, G.W.: An improved elastic stain. Bull. Inst. Med. Lab. Technol., 9, 9, 1943
Van Gieson, I.: Laboratory notes of technical methods for nervous system. New York, Med. J., 50:57- 60; 1889
Verrhoeff, F. H.: Some new staining methods of wide applicability. Including a rapid differential stain for elastic tissue. J. Amer. Med. Ass., 50, 876, 1908
Weigert, C.: Ueber einer Methode zur Färbung elastischer Fasern. Zbl. allg. Path. Anat., 9, 289, 1898

F01IT04_V1R0412_AP0384_English



Catalog number



Batch code



In Vitro diagnostic medical device



Temperature limitation



Expiration date



Test number



Manufacturer



See instruction for use



Genova Scientific, S.L.
C/ Johann Gutenberg, 4F. Pol. Ind.
El Cádizamo I • 41300 San José
de La Rinconada • Sevilla, SPAIN
Teléfono: +34 954 150767
Fax: +34 955 266494

info@gennovalab.com
www.gennova-europe.com