

EcoNova Taq DNA Polymerase

Without dNTPs

Reference: AB12000; AB12001; AB12002



1 of 2

INTENDED USE AND PRESENTATION:

This polymerase is recommended for use in molecular biology.

AB12000, 500 units. Concentration of enzyme 5 units/ μ L.

AB12001, 1000 units. Concentration of enzyme 5 units/ μ L.

AB12002, 2000 units. Concentration of enzyme 5 units/ μ L.

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 η moles of dNTPs into acid-insoluble material in 30 minutes at 74°C.

For research use only.

SUMMARY, EXPLANATION AND LIMITATIONS:

EcoNova Taq DNA Polymerase is a thermostable high quality recombinant enzyme derived from an *E. coli* strain that carries an overproducing plasmid containing a modified gene of *Thermus aquaticus* DNA Polymerase. It is able to withstand repeated heating to 95°C without significant loss of activity and the estimated half life at 95°C is 1,5 hours. The enzyme is ~94 kDa by SDS-PAGE catalyzing 5'→3' synthesis of DNA. The enzyme has no detectable 3'→5' proofreading exonuclease activity. It is the latest developments in polymerase technology and buffer chemistry to enhance PCR speed, yield and specificity. The enzyme and buffer system allow for superior PCR performance on complex templates such as mammalian genomic DNA.

EcoNova Taq DNA Polymerase is provided with 10 fold concentrated (X) reaction buffer that contains PCR enhancers. This reaction buffer will enable or improve sub-optimal PCR caused by templates that have a high degree of secondary structure or that are GC-rich.

EcoNova Taq DNA Polymerase has an error rate of approximately 1 error per 2,5x10⁵ nucleotides incorporated.

This product is sold for research purposes only. All claims must be brought within expired date.

Quality Control: Endonuclease, Exonuclease, DNase, RNase and Protease activity is not detected.

EcoNova Taq DNA Polymerase is determined to be >98% pure as judged by SDS-PAGE.

APPLICATIONS:

The EcoNova Taq DNA Polymerase is a robust enzyme for all your everyday PCR applications including genotyping, screening and library construction. This Taq DNA Polymerase can perform consistently well on a broad range of templates (including both GC and AT rich).

PRODUCT COMPOSITION:

-10X Reaction Buffer: Contains 800 mM Tris-HCl (pH 9,0), PCR enhancers, (NH₄)₂SO₄, 0,2% w/v Tween-20.

-25 mM MgCl₂.

METHODS AND PROCEDURE:

Optimal reaction conditions, such as reaction time, temperature, and amount of template DNA, may vary and must be individually determined.

General Reaction Protocol:

1. Thaw 10 fold concentrated reaction buffer and dNTPs mixture.
2. Prepare a master mix.

Component	Volume	Final Conc.
10X Reaction Buffer	10 μ L	1X
20 mM dNTPs Mixture	1 μ L	200 μ M
25 mM MgCl ₂	6~10 μ L	1,5~2,5 mM
Upstream Primer	1~3 μ L	0,1~0,3 μ M
Downstream Primer	1~3 μ L	0,1~0,3 μ M
EcoNova Taq DNA Polymerase (5 U/ μ L)	0,4~1,0 μ L	0,02~0,05 U/ μ L
Template DNA	5~20 μ L	5~100 η g/ μ L
Sterilized D.W.	Up to 100 μ L	-
Total Volume	100 μ L	-

*Primers should have a predicted melting temperature of around 60°C, using default Primer 3 settings (<http://frodo.wi.mit.edu/primer3/>).

Amount of template:

Bacteriophage λ , cosmid, plasmid → 10 fg~300 η g.

Total genomic DNA → 100 η g~1 μ g.

3. Mix the master mix and dispense appropriate volumes into PCR tubes. Centrifuge the reactions in a microcentrifuge for 10 seconds.
4. Perform PCR using your standard parameters (3-step).

Step	Temperature & Reaction Time		
Initial denaturation	3~5 min. at 94°C	-	-
25~35 cycles	30~60 sec. at 95°C	30~60 sec. at 50~68°C	1~4* min. at 72°C
Final extension	-	-	5~10 min at 72°C

*For PCR products longer than 3~4 Kb, use an extension or elongation time of approximately 1 min per Kb DNA.

5. Separate the PCR products by agarose gel electrophoresis and visualize with EtBr or any other means.

A DNA fragment which is amplified by Nova Taq DNA Polymerase has A-overhang, and it enables you to do cloning by using T-vector.

REQUIRED MATERIALS BUT NOT SUPPLIED:

All reagents, materials, and laboratory equipment for PCR procedures are not provided with this polymerase. This includes sterile reaction tubes, micropipettes and tips, template DNA, gen-specific PCR primer pair, PCR grade H₂O, heat pretreatment equipment (thermoblock, microwave), centrifuge, cold store and thermal block cycler.

Buffered solutions for DNA extraction or purification, enzyme treatments, and other auxiliary reagents are available from Genova Scientific.



Catalog number



Batch code



Research use only



Temperature limitation



Expiration date



Manufacturer



See instruction for use



Genova Scientific, S.L.
C/ Johann Gutenberg, 4F. Pol. Ind.
El Cafamo I • 41300 San Jose
de La Rinconada • Sevilla, SPAIN
Telefono: +34 954 150767
Fax: +34 955 266494

info@gennovalab.com
www.gennova-europe.com

EcoNova Taq DNA Polymerase

Without dNTPs

Reference: AB12000; AB12001; AB12002



2 of 2

STORAGE AND STABILITY:

Store at -20°C until the expiration date printed on product label. Shipping and temporary storage for up to 1 month at room temperature has no detrimental effects on the quality of EcoNova Taq DNA Polymerase. Avoid prolonged exposure to light. Avoid multiple freeze-thaw cycles and exposure to frequent temperature changes. Do not use it after the expiration date. If the product is stored under different conditions from those stipulated in these technical indications, the new conditions must be verified by the user. The validity period of the ready to use products when opened, is the same as the expiration date indicated on the label of intact product.

Gennova Scientific guarantees that the product will maintain all of the described characteristics from the production date until the expiration date, as long as the product is stored and used as recommended. No other guarantees are provided. Under no circumstances Gennova Scientific is obliged to cover damages caused by use of this reagent.

TROUBLESHOOTING:

If unusual amplification is observed or any other deviations from the expected results, please read these instructions carefully, along with the instructions from the PCR system. If this does not solve the problem, please contact Gennova Scientific's technical support department (techsupport@gennovalab.com) or your local distributor.

PRECAUTIONS:

Use proper protective equipment in order to avoid contact with reagents and samples in the eyes, skin, and mucosal tissues. In case of contact with sensitive areas, immediately flush the affected area with water. Avoid microbial contamination of the reagent, as this may produce nonspecific amplification results.

Use only by qualified personnel.

Material safety data sheet (MSDS) is available upon request.

PERFORMANCE CHARACTERISTICS:

Gennova Scientific has performed studies to evaluate the functioning of this polymerase for use with standard amplification systems, concluding that the product is both specific and sensitive for PCR performance.

BIBLIOGRAPHY:

Chien A., Edgar D.B., Trela J.M., "Deoxyribonucleic acid polymerase from the extreme thermophile *Thermus aquaticus*", *Journal of Bacteriology*, 127(3), 1550-57, 1976.
Lawyer F.C., Stoffel S., Saiki R.K., Myambo K., Drummond R., et al., "Isolation, characterization, and expression in *Escherichia coli* of the DNA polymerase gene from *Thermus aquaticus*", *The Journal of Biological Chemistry*, 264(11), 6427-37, 1989.
Tindall K.R., Kunkel T.A., "Fidelity of DNA synthesis by the *Thermus aquaticus* DNA polymerase", *Biochemistry*, 27(16), 6008-13, 1988.
Innis M.A., Myambo K.B., Gelfand D.H., Brow M.A., "DNA sequencing with *Thermus aquaticus* DNA polymerase and direct sequencing of polymerase chain reaction-amplified DNA", *Proceedings of the National Academy of Sciences of the United States of America*, 85(24), 9436-40, 1988.
Lo Y.M., Mehal W.Z., Fleming K.A., "Rapid production of vector-free biotinylated probes using the polymerase chain reaction", *Nucleic Acids Research*, 16(17), 8719, 1988.
Erich H.A., (ed.) 1988, "PCR technology: principles and applications for DNA amplification", Stockton Press, New York.

F01IT04_AB12000_AB12001_AB12002_V1R0912_EN_EcoNova_Taq_DNA_Polymerase



Catalog number



Batch code



Research use only



Temperature limitation



Expiration date



Manufacturer



See instruction for use



Gennova Scientific, S.L.
C/ Johann Gutenberg, 4F. Pol. Ind.
El Cafamo I • 41300 San Jose
de La Rinconada • Sevilla, SPAIN
Telefono: +34 954 150767
Fax: +34 955 266494

info@gennovalab.com
www.gennova-europe.com