

Anti-IDH1 R132H / DIA-H09

Mouse monoclonal anti-brain tumor marker (Astrocytoma, Oligodendroglioma), Clone H09

Product Information

Catalog No.:	DIA-H09 (100µg) DIA-H09-M (20µg)	Reconstitution:	DIA-H09 (100µg), restore to 500 µl DIA-H09-M (20µg), restore to 100 µl Reconstitute with sterile distilled water by gentle shaking for 10 minutes
Clone:	H09	Presentation:	In PBS with 2% BSA, 0.05% NaN ₃ , pH 7.4. Antibody purified from culture supernatant
Concentration:	0.2 mg/ml	Applications:	Immunohistochemistry (standard formalin-fixed paraffin sections) Western blot
Isotype:	Mouse IgG2a	Dilutions:	1:20 Immunohistochemistry (IHC) 1:500 Western Blot (General recommendation, validation of antibody performance/protocol is the responsibility of the end user. Positive/negative controls should be run simultaneously with patient specimen. Interpretation must be made by a qualified pathologist within the context of patient's clinical history/other diagnostic tests.)
Specificity:	Human IDH1 R132H point mutation	Control Antibody:	DIA-W09, rat anti-hu IDH1wt, clone W09
Immunogen:	Synthetic peptide, amino acid sequence CKPIIIGHHAYGD		
Physical State:	Lyophilized powder		
Species	Human		
Reactivity:	Human		
Positive Control:	Oligodendroglioma, diffuse astrocytoma		
Negative Control:	Pilocytic astrocytoma, primary glioblastoma (ca. 95% of cases negative)		
Visualization:	Cytoplasmic		

(Staining protocols on www.dianova.com)

Reactivity

Antibody clone H09 reacts specifically with the isocitrate dehydrogenase 1 (IDH1) R132H point mutation in tissue sections from formalin-fixed brain tumor specimens. Heterozygous point mutations of IDH1 codon 132 are frequent in World Health Organization (WHO) grade II and III gliomas. IDH1 R132H mutations occur in approximately 70% of astrocytomas and oligodendroglial tumors. The high frequency and distribution of the IDH1 R132H mutation among specific brain tumor entities allow the highly sensitive and specific discrimination of various tumors by immunohistochemistry, such as anaplastic astrocytoma from primary glioblastoma or diffuse astrocytoma WHO grade II from pilocytic astrocytoma or ependymoma. Noteworthy is the discrimination of the infiltrating edge of tumors with IDH1 mutation from reactive gliosis. This antibody is highly useful for tumor classification and in detecting single infiltrating tumor cells.

Instructions for Use

Immunohistochemical staining of standard formalin-fixed paraffin sections

Deparaffinize and rehydrate according to standard procedures. Heat induced epitope retrieval (HIER) is required. For immunohistochemical detection different techniques can be used: Indirect immunoenzyme labeling with a secondary antibody conjugate, biotin(strept)avidin-based detection, soluble enzyme immune complex or polymer-based detection. To detect antibody, follow the instructions provided with the particular visualization system. The antibody is suited for immuno-histochemical staining using automated platforms. Use the antibody at 1:20 dilution for 30min at RT.

Technical note

Diffuse astrocytoma WHO grade II may have low protein-expression. At high dilution of the antibody single tumor cells in the infiltration zone may not be stained.

Intended use / regulatory status

Europe: For in Vitro Diagnostic Use / All other countries: For Research Use only

Storage and Stability

The H09 antibody in lyophilised form is stable for at least one year. As reconstituted liquid store at 2-8°C short term (several weeks). For long term storage aliquot and freeze at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Safety Notes

The material contains 0.05% sodium azide as preservative. Although the quantity of azide is very small, appropriate care should be taken when handling this material. Avoid skin and eye contact, inhalation, and ingestion.




Figures

Immunohistochemistry of human IDH1 R132H in formalin-fixed paraffin-embedded brain tissue sections

(pictures courtesy of Prof. Dr. med. Andreas von Deimling, Department of Neuropathology, University Heidelberg / Clinical Co-operation Unit Neuropathology, German Cancer Research Center (DKFZ), Heidelberg, Germany)

- A:** Strong reaction of IDH1 mutation specific antibody clone H09 in tumor center of anaplastic oligoastrocytoma.
- B:** Infiltration zone of anaplastic astrocytoma with specific labelling of infiltrating glioma cells by antibody clone H09.
- C:** Identification of single tumor cells in white matter distant from tumor center with IDH1 mutation specific antibody clone H09.
- D:** Cortex infiltrated by oligodendroglioma with specific labelling of tumor cells by antibody clone H09.
- E:** Double staining of GFAP (glial fibrillary acidic protein, red) and clone H09 (brown) of oligodendroglioma infiltration zone demonstrating specific labelling of tumor cells but not of GFAP positive reactive astrocytes.
- F:** Strong reaction of IDH1 mutation specific antibody clone H09 with IDH1 R132H mutated diffuse astrocytoma (left) but not with wild type tumor (right).



References

1. Capper D, Zentgraf H, Balss J, Hartmann C, von Deimling A. Monoclonal antibody specific for IDH1 H132R mutation. *Acta Neuropathol.* 118(5): 599-601, 2009
2. Capper D, Weissert S, Balss J, Habel A, Meyer J, Jäger D, Ackermann U, Tessmer C, Korshunov A, Zentgraf H, Hartmann C, von Deimling A. Characterization of R132H mutation-specific IDH1 antibody binding in brain tumors. *Brain Pathol.* 20(1): 245-254, 2010
3. Capper D, Sahm F, Hartmann C, Meyermann R, von Deimling A, Schittenhelm J. Application of mutant IDH1 anti-body to differentiate diffuse glioma from nonneoplastic central nervous system lesions and therapy-induced changes. *Am J Surg Pathol.* 34(8):1199-1204, 2010

Symbols

REF	Catalog Number		Expiry	CE	Conformity with IVDD 98/79/EC
LOT	Lot Number		Temperature Limitation	IVD	For In vitro Diagnostic Use
	Manufacturer		Consult Instructions for Use		